

Traceability Statement for Urine analyzer Calibration liquid for Specific gravity

Product name: Urine analyzer Calibration liquid for Specific gravity

1. Traceability Description:

In accordance with ISO 17511 and CLSI regulations, the value of urine analyzer Calibration liquid for Specific gravity is strictly traced back to the refractive specific gravity meter.

2. Value Traceability Procedure:

According to ISO 17511 "In vitro diagnostic medical devices - Measurement of quantities in biological samples - Metrological traceability of values assigned to calibrators and control materials," the value of urine analyzer Calibration liquid for Specific gravity is strictly traced back to the refractive specific gravity meter.

