

iQM2 Corrective Action Report

November 2025

Analyzer

Model: GEM® Premier 5000
Area: POCT15TALI
Name: GEM-5000
S/N: 20075305

| Cartridge | Date/Time | Event | Corrective Action | Result |
|--|--|--|--|------------------------|
| Lot No.: 250716H S/N: 0000000005012 37510 Exp. Date: 19.12.2025 | 24.11.2025 02:49:59 18.11.2025 11:08:43 | Solution A Error for Glu Solution A Error for Cl ⁻ | Sensor Output Adjusted Sensor Output Adjusted | Corrected Corrected |
| Lot No.: 250716H S/N: 0000000005012 37507 Exp. Date: 18.11.2025 | 18.11.2025 07:55:45 04.11.2025 00:32:17 | Cartridge Removed. Samples Remaining = 14. No. of Solution B Adjustments = 74 Solution D Error for O ₂ Hb, COHb, MetHb, HHb, tHb | Sensor Output Adjusted | Corrected |