

iQM2 Corrective Action Report

November 2025

Analyzer

Model: GEM® Premier 5000
Area: POCT06KCLI
Name: GEM5
S/N: 19042829

Cartridge	Date/Time	Event	Corrective Action	Result
Lot No.: 250716H S/N: 0000000005012 37492 Exp. Date: 08.12.2025	26.11.2025 02:51:30	Solution A Error for pO_2	Sensor Output Adjusted	Corrected
	24.11.2025 02:50:01	Solution A Error for Glu	Sensor Output Adjusted	Corrected
	17.11.2025 02:50:03	Solution A Error for Glu	Sensor Output Adjusted	Corrected
Lot No.: 250515W S/N: 0000000005011 81685 Exp. Date: 09.11.2025	07.11.2025 11:58:31	Cartridge Removed. Samples Remaining = 10. No. of Solution B Adjustments = 26		