

Supplier:
 Terumo BCT, Inc.
 10811 West Collins Avenue
 Lakewood, CO 80215, USA
 Toll Free (US & Canada): 877-339-4228
 Phone :303-231-4357
 FAX :303-542-5215

Manufacturer:
 TERUMO CORPORATION Fujinomiya Factory
 818, Misonodaira, Fujinomiya-shi, Shizuoka 418-0004 JAPAN
 Phone :+81-544-27-7142
 FAX :+81-544-23-4548
 Quality Assurance Department:
 Yoshito Ide

LOT CERTIFICATE


Product Name:	TERUFLEX TRANSFER BAG		
Product Code/Cat No:	BB*T100BM		
Lot number:	250721F3		
Lot size (amount):	9,600 sets		
Sterilization Date:	22/07/2025	Sterilization Method:	Ethylene Oxide
Expiry:	06/2028		

Material validation Quality of product is validated according to ISO3826-1.

QC specifications and inspection result	(1)Chemical test for PVC material based on internal standards	<input checked="" type="checkbox"/>	Passed
	(2)Chemical test for other plastic material based on internal standards	<input checked="" type="checkbox"/>	Passed
	(3)Function test for PVC container based on internal standards	<input checked="" type="checkbox"/>	Passed
	(4)Biological test (Bacterial Endotoxin) according to ISO 3826-1	<input checked="" type="checkbox"/>	Passed
	(5)Microbiological test (Biological Indicator test) according to ISO 11135 The sterilization assurance level of a validated sterilization process is 1 or less visible microorganisms per 1x10E06 finished products	<input checked="" type="checkbox"/>	Passed

Comments results: Lot No. 250721F3 meets all requirements of QC product
 Product Code/Cat No: BB*T100BM specification

Release statement : Lot Released
 Date of release determination : 08 / 08 / 2025
 (Day) (Month) (Year)



 Yoshito Ide,
 Quality Assurance Department
 Terumo Corporation Fujinomiya Factory

Date of signature : 21 / 10 / 2025
 (Day) (Month) (Year)